濫用藥物確認檢驗判定標準 (尿液) By Shih


Drug Class	Analyte	Minirtry of Health and Welfare Cutoff (Taiwan)	The Substance Abuse and Mental Health Services Administration Cutoff (USA)
Amphetamines 安非他命類藥物			
Amphetamines 安非他命	Amphetamines	500 ng/mL	250 ng/mL
Methamphetamine 甲基安非他命	Methamphetamine	500 ng/mL	250 ng/mL
	Amphetamines	100 ng/mL	
Methylenedioxymethamphetamine (MDMA) 3,4-亞甲基雙氧甲基安非他命	MDMA	500 ng/mL	250 ng/mL
	MDMA+MDA	500 ng/mL	
Methylenedioxyamphetamine (MDA) 3,4-亞甲基雙氧安非他命	MDA	500 ng/mL	250 ng/mL
Methylenedioxyethylamphetamine (MDEA) 3,4-亞甲基雙氧-N- 乙基安非他命	MDEA	500 ng/mL	250 ng/mL
Opiate Metabolites 海洛因、鴉片代謝物	Codeine	300 ng/mL	2000 ng/mL
	Morphine	300 ng/mL	2000 ng/mL
Marijuana Metabolites 大麻代謝物	Delta-9-tetrahy drocannabinol-9 carboxylicacid (THCA)	15 ng/mL	15 ng/mL
Cocaine metabolites 古柯鹼代謝物	Benzoylecgonine	150 ng/mL	100 ng/mL
Ketamine metabolites 愷他命代謝物	Ketamine	100 ng/mL	-
	Norketamine 去甲基愷他命	100 ng/mL	_
	Ketamine+Norketamine	100 ng/mL	_


Reference: 濫用藥物尿液檢驗作業準則,政院衛生署署授管字第0960710315號; Federal Register, November 25, 2008 (73 FR 71858), Section 3.4



Dilute and Shoot

SQUARESPACE.COM/LOGO ICONS BY THE NOUN PROJECT





Extraction or HPLC

Strata-X-Drug B is a polymeric strong cation-exchange SPE sorbent that is designed and quality controlled with Substance Abuse and Mental Health Services Administration (SAMHSA) regulations in mind.

Strata-X-Drug B also Provides:

- Time and solvent savings; no conditioning required
- Consistent, high recoveries of SAMHSA regulated drugs
- Sensitivity at SAMHSA's newly regulated cutoff levels

